
# Panasonic®

## Operating Instructions 使用說明書

### (Household) Air Massager (家用) 空氣按摩器

Model No.    EW-RAJ6  
型號

Leg Refre

English	EN2
繁體中文	TC2

Panasonic Corporation  
© Panasonic Corporation 2017  
<http://www.panasonic.com>

No.2 EN, TC

Printed in China  
中國印刷  
W9030RAJ422 F1012Y1

Thank you for purchasing this Panasonic product.

**Before operating this unit, please read these instructions completely and save them for future use.**

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# What is the Air Massager?

The Air Massager is the household massager that uses airbags for pressure massage. Wrap your legs with the attachments and turn on the manipulator. The airbags swell and massage your legs.


In this case,  
use for 10 minutes  
a day

Program	For the following cases
<b>Toe ~ Knee</b>	
	<div style="display: flex; justify-content: space-between;"> <div style="width: 25%;">All in one</div> <div>When you want to refresh your toes to around knees.</div> </div> <div style="display: flex; justify-content: space-between;"> <div>Squeeze</div> <div>When your legs feel heavy.</div> </div> <div style="display: flex; justify-content: space-between;"> <div>Knead*</div> <div>When you want to relieve your muscle tension by comfortable massage.</div> </div> <div style="display: flex; justify-content: space-between;"> <div>Toe ~ Calf</div> <div>When your toes and calves are swelling.</div> </div> <div style="display: flex; justify-content: space-between;"> <div>Rear Knee</div> <div>When you cannot recover from the whole leg fatigue and tiredness.</div> </div>
<b>Toe ~ Thigh</b>	
	<div style="display: flex; justify-content: space-between;"> <div style="width: 25%;">All in one</div> <div>When you want to refresh your toes to thighs.</div> </div> <div style="display: flex; justify-content: space-between;"> <div>Soft Knead*</div> <div>For relaxing before going to bed or muscle strain after exercise.</div> </div> <div style="display: flex; justify-content: space-between;"> <div>Deep Tissue*</div> <div>When you feel tiredness in your thighs.</div> </div>

• Refer to Page EN20-EN21 for details of the program and massage.

\* Pressure function like knead

**Warmth**

**Warmer (heater)**

- The warmer (heater) warms the textile around toes only.

## Contents

**Safety Precautions** EN4

**Name of each component** EN8

**Selected part and methods of usage** EN10

**Preparations** EN14

- Wrap it around the knee
- Wrap it around the thigh
- Use the sole pads / shin pads

**Massage** EN20

- To change the program
- Adjust temperature around toes
- To change the intensity

**After use** EN22

- Storage
- Keeping

**Cleaning and maintenance** EN23

**Frequently asked questions** EN24

**Before entrusted repair** EN26

**Specifications** EN28

# Safety Precautions

Always follow these instructions

This section contains instructions that you must follow to prevent personal injury or damage to property.

■ **The instructions are classified according to the level of injury or damage that may occur due to the mistaken use of this unit. Please ensure that you read all instructions before using the Air Massager.**

**⚠ WARNING** Denotes a potential hazard that could result in serious injury or death.

**⚠ CAUTION** Denotes a hazard that could result in minor injury or property damage.

## ⚠ WARNING

### Symptoms

● **The following persons should not use this unit.**

- (1) Persons who are prohibited from receiving massages by a physician (e.g.: due to thrombosis, severe aneurism, acute varicose veins, any type of dermatitis and skin infections (including inflammation of the hypodermis), etc.)
- (2) Persons suffering from osteoporosis, spine fracture, sprain or acute pain such as a pulled or torn muscle
- (3) Persons who have a pacemaker, or other personal medical devices prone to electronic interference
- (4) Pregnant women and women who have just given birth
- (5) Persons suffering from severe circulatory disorders in the legs (Otherwise it may cause an accident, ill health or worsening of symptoms.)

● **This unit is not intended for use by persons (including children) with reduced physical sensory or mental capabilities, or lack of experience and knowledge, unless they have been given supervision or instruction concerning use of the unit by a person responsible for their safety. Children should be supervised to ensure that they do not play with the unit.**

(Doing so may cause an accident or injury.)

● **Persons with any of the conditions listed below or who are undergoing medical treatment should consult with their physicians before using the unit.**

- (1) Persons who have a malignant tumor
- (2) Persons suffering from heart disease
- (3) Persons who have a loss of heat sensation
- (4) Persons who have sensory abnormalities due to a severe peripheral circulatory disturbance such as diabetes mellitus
- (5) Persons who have a wound where the massager operates
- (6) Persons whose body temperature is over 38 °C (who have a fever) (E.g.: When a person has acute inflammation symptoms (fatigued, chills, blood pressure fluctuations, etc.) or when debilitated)
- (7) Persons who require bed rest or who are in poor physical shape
- (8) Persons other than those listed above who feel unwell

● **If you begin to feel unwell during use or if you do not feel effects of the air massager, stop use immediately and consult a physician.**

■ **Instructions that you must follow are indicated using the following symbols.**  
(The following symbols are examples.)

Denotes a specific operating procedure that must not be performed.

Denotes a specific operating procedure that must be performed.

● **The unit has a heated surface. Persons insensitive to heat must be careful when using the unit.**

(Failure to do so may cause an accident, injury or become unwell.)

### When in use

● **Do not use the right attachment for left, and vice versa.**  
(Otherwise it may lead to poor body conditions.)

● **Do not use sole or shin pads for parts other than those specified or do not apply them to the bone.**  
(Otherwise it may lead to injuries.)

● **Do not stand or walk when in use.**

● **Do not use the unit only for one leg.**

● **Don't use for purposes other than massaging around the feet, calves, around knees and thighs.**

● **Do not allow children to use the unit, make sure they do not play with the manipulator or attachment.** (Doing so may cause an accident or injury.)

● **Do not use the unit in bedding or kotatsu.**

(Otherwise it may lead to burn or electric shock.)

● **The service time shall be less than 10 minutes per day.**

(Otherwise it may lead to a low-temperature burn\*, the opposite effect or injuries.)

\* Burns may occur at relatively low temperatures (40 °C to 60 °C) when in contact with the same area of skin for long periods of time, even if you do not feel heat or pain.

The elderly, persons having a delicate skin or insensitive cutaneous sensation may have redness on their skin just after use. If the redness does not disappear after one hour, consult the doctor.

● **Before use, make sure that the cloth is not damaged. Stop using, unplug the AC adaptor and entrust the repairs to qualified persons even for minimum damages.**

(Otherwise it may lead to injuries or electric shock.)

● **For initial use, select "Intensity: Low" mode until you are accustomed to it.** (Otherwise it may lead to injuries or poor body conditions.)

# Safety Precautions (Continued)

## ⚠ WARNING

### Power plug, power cord etc.

- **Do not do anything that may damage the power cord or power plug.**  
Do not scratch, modify, use near a heater, excessively bend, twist, pull, place heavy objects on top of, or bundle the cord or plug.  
(Doing so may cause electric shock, burning or fire due to a short circuit.)  
→ Contact the nearest authorized service center for repairs to the power cord or power plug.
- **When you unplug the AC adaptor, do not pull the power cord.**  
(Doing so may cause electric shock, burning or fire due to a short circuit.)
- **Do not plug or unplug the AC adaptor or air tube plug with wet hands.**  
(Doing so may cause electric shock.)

- **Always use sockets and wiring devices at the correct rated value.**
- **Use 220 V – 240 V ~ (Cannot be used overseas, nor with a transformer.)**
- **Always insert the power plug fully into the socket.**  
(Failure to do so may result in electric shock or fire due to overheating.)

### If a malfunction or breakdown occurs


- **If a malfunction or breakdown occurs, stop using the unit immediately and unplug the power plug from the wall socket.**  
(Failure to do so may result in smoke / fire or electric shock.)

<Examples of malfunctions / breakdowns>

- AC adaptor or air tube plug is damaged
- The plug of the AC adaptor is loose
- The air tube plug is loose when it is inserted to the manipulator

→ Please do not repair it yourself. Entrust the authorized service center with inspection and repairs immediately.

### Take note of the following points as well

- **Do not stick a pin or needle into the attachment or scar with a knife.**  
(Otherwise it may cause fire, an electric shock or burn due to damaged warmer (heater) on the attachment.)
- **Do not modify or disassemble.** (Doing so may cause fire, or may cause the unit to function abnormally resulting in injury.)

## ⚠ CAUTION

### Symptoms

- **If symptoms such as rashes, reddening, itching, etc. occur as a result of using this unit, stop using the unit and consult a physician.**  
(Failure to do so may cause an accident or adversely affect health.)

### Power supply


- **During cleaning and maintenance, unplug the AC adaptor from the socket.** (Failure to do so may cause an electric shock or burn.)
- **Unplug the AC adaptor from the socket when not in use.** (Since dust and humidity deteriorates insulation, which may result in fire due to electrical fault.)
- **The unit is only to be used with the AC adaptor provided with unit.**

### Before and during use

- **Don't use in places with high humidity, such as bathrooms.**  
(Doing so may cause electric shock.)
- **Do not immerse manipulator or attachment in water or wash with water.** (Doing so may cause electric shock, burning or fire due to a short circuit.)
- **Do not fall asleep while using the unit.**
- **Do not use the unit after drinking alcohol.** (Doing so may cause an accident or injury.)
- **Don't allow pins, garbage or water etc. to enter the AC adaptor's power plug or the gap of the appliance plug.**
- **When using the massager, don't unplug the AC adaptor from the socket or pull out the appliance plug from the manipulator side.**  
(Doing so may cause electric shock, burning or fire due to a short circuit.)
- **When pulling out the air tube plug, do not pull the air tube for connection.**
- **For storage, do not wrap the power cord around the manipulator or the power plug.**
- **Do not fold the air tube for connection or wrap it on the attachment.** (Otherwise, the power cord or air tube for connection may be disconnected due to an overload of the power cord, causing fire or electric shock due to short circuit.)
- **Do not use at the same time as another medical device.**  
(Doing so may cause an accident or deteriorated health.)
- **Don't use accessories other than the provided sole and shin pads.** (Otherwise it may lead to injuries.)
- **Do not use AC adaptors other than those specified (Model No.: EW-3A08) or use the provided adaptor for other equipment.**  
(Doing so may cause electric shock, burning or fire due to a short circuit.)


- **Unplug the AC adaptor immediately in case of a power failure.**  
(Otherwise it may lead to accidents.)

# Name of each component

## Manipulator


The following symbol indicates that a specific detachable AC adaptor is required for connecting the electrical appliance to the supply mains. The type reference of AC adaptor is marked near the symbol.


## Accessories

### Sole pads (2 pcs)

Protrusion

### Shin pads (2 pcs)


Handlebar

Protrusion

## Attachment


The unit is composed of right and left attachments. (Bilateral symmetry)  
The illustration shows the left attachment.



Name of each component

## Selected part and methods of usage

Airbags arranged to wrap the whole leg provides comfortable massage for the target point.


You can enjoy a comfortable massage with **12 massaging points ★** over the whole leg with two wrapping varieties, "Toe ~ Knee" and "Toe ~ Thigh".

## Selected part and methods of usage (Continued)


You can enjoy two wrapping varieties, "Toe ~ Knee" and "Toe ~ Thigh".

### Toe ~ Knee

#### Above knee (★7)

Kneading massage like a strong-grip hand massage.

#### Below knee (★9)

Loosen the stiffness of the whole leg by loosening the joints of calf and thigh muscles.


#### Above calf (★8)

Massage the upper calf like it was being gripped.

Apply the point ball to the back of the knee for a compression massage

Reproduces the comfort of pressing from a human palm. Loosen the back of the knees with the airbags.

### Toe ~ Thigh

#### Front thigh (★10)

Firmly loosen the quadriceps femoris muscle that is mainly on the side of the thigh.

#### Backside of the thigh (★12)

Massage large hamstring muscle on the backside of the thigh firmly.

#### Inside / outside of the thigh (★11)

Strong massage, as if the inside / outside of the thigh was being gripped.

The point ball firmly pushes the thigh

Loosen the hamstring muscle on the backside of the thigh firmly like pressing with the palm.

# Preparations

## Wrap it around the knee

### First of all, confirm the following contents

- Whether the cloth is damaged?
- Whether it operates normally after switching on? (The confirmation is also necessary if it is not used for a certain period of time)
- In case of abnormalities, please refer to Page EN26-EN27.

#### • The unit may leave marks after use. If you are concerned about marks...

Put the unit on over socks and trousers.

The unit leaves marks on calves / feet / around knees and they may not disappear immediately. If you use the unit before going out, note this point.

#### • Be careful not to have the surface adhesive tape contact with your clothes or stockings when putting on the unit.

#### • If you put on the unit over thick trousers or leg warmers, the surface adhesive tape may be undone.

#### • If you want a stronger massage...


We recommend that you use sole and shin pads (Page EN18-EN19)

#### • If you want to feel stronger compression from the point balls...

Try it with your legs stretched. (Page EN21)


## 1 Connect the AC adaptor

220 V – 240 V ~

## 2 Check the right and left attachments and the position of point balls before you wrap the attachments



- 1) Check the right and left attachment by the label
- 2) Check the position of the point balls by touching them with a hand

## 3 Line your heels into the internal Receiver (in the textile) and attach it in sequence from the feet.

e.g.) for left

- 1) Hold the end of the surface adhesive tape and fix it in order of your **Toe** → **Ankle** → **Calf**
- 2) Press the end of the surface adhesive tape again firmly so that it does not loosen

## 4 Wrap it around the knee


- 1) Touch with your hands, check the point ball and align to the back of your knees  
You can overwrap the strap
- 2) Wrap it to squeeze your knees  
• Fix the straps in order of **Below knee** → **Above knee**
- 3) Wrap the right attachment as in procedures 3 to 4

• You can adjust the intensity of massage by how you wrap the attachments around your legs. (Page EN25)

## 5 Insert the air tube plug into the manipulator

After the preparations are completed...

Massage (Page EN20)

Preparations (Wrap it around the knee)

# Preparations

## Wrap it around the thigh

First of all, confirm the following contents


- Whether the cloth is damaged?
- Whether it operates normally after switching on? (The confirmation is also necessary if it is not used for a certain period of time)
- In case of abnormalities, please refer to Page EN26-EN27.

- **The unit may leave marks after use. If you are concerned about marks...**  
Put the unit on over socks and trousers.  
The unit leaves marks on calves / feet / around knees / thighs and they may not disappear immediately. If you use the unit before going out, note this point.
- **Be careful not to have the surface adhesive tape contact with your clothes or stockings when putting on the unit.**
- **If you put on the unit over thick trousers or leg warmers, the surface adhesive tape may be undone.**
- **If you want a stronger massage...**  
We recommend that you use sole and shin pads (Page EN18-EN19)


### 1 Connect the AC adaptor

220 V – 240 V ~


### 2 Check the right and left attachments by the label before you wrap the attachments


### 3 Line your heels into the internal Receiver (in the textile) and attach it in sequence from the feet.

e.g.) for left

- 1) Hold the end of the surface adhesive tape and fix it in order of your **Toe** → **Ankle** → **Calf**
- 2) Press the end of the surface adhesive tape again firmly so that it does not loosen


Receiver

Preparations (Wrap it around the thigh)

### 4 Wrap it around the thigh

- 1) Wrap the leg with the strap above the knee
- 2) Fix the straps in order of **Thigh** → **Above knee**
  - You can easily wrap the strap while sitting on the edge of a chair.


3) Wrap the right attachment as in procedures 3 to 4

- You can adjust the intensity of massage by how you wrap the attachments around your legs. (Page EN25)

### 5 Insert the air tube plug into the manipulator

After the preparations are completed...

Massage (Page EN20)

# Preparations



## Use the sole pads / shin pads

**If you want a stronger massage** When you are tired after exercise or walking with high-heeled shoes...

### Sole pads

**1** Preparation before setting  
(Procedure 1 to 2 on Page EN14, 16)


**2** Set the sole pads to the pockets

•Your may have a pain in your hand if it touches the surface adhesive tape.


#### Plantar arch

Adjust the large protrusion inside to suit your sole

#### Other (comfortable parts)

Locate a comfortable part and adjust the protrusion to it

**3** Wrap the attachments around your legs  
(Procedure 3 to 4 on Page EN15, 17)

**4** Insert the air tube plug into the manipulator  
(Procedure 5 on Page EN15, 17)

- When you use the sole pads, you may be less likely to feel warmth.
- The height of the protrusions are represented with large and small.


### Shin pads

**1** Preparation before setting  
(Procedure 1 to 2 on Page EN14, 16)

**2** Locate the most comfortable position after pressing


- Press around the shin muscles with the shin pad and locate the most comfortable position.

**3** Wrap the attachment around your leg from the toe to calf  
(Procedure 3 on Page EN15, 17)

**4** Insert the shin pad straightly into the part you locate in procedure 2

- Insert the protrusion toward the muscles
- Be careful not to touch the bone
- Set the right and left shin pads in the same way

**5** Wrap the attachment straps around your knee or thigh  
(Procedure 4 on Page EN15, 17)

**6** Insert the air tube plug into the manipulator  
(Procedure 5 on Page EN15, 17)

Preparations (Use the sole pads / shin pads)

# Massage

•To change the program

•Adjust temperature around toes

•To change the intensity

Standard massaging time The service time shall be less than 10 minutes per day

- After using for a long time, it will produce excessive stimulation against the muscles or nerves and sometimes it will lead to the opposite effect.

## By pressing , it will switch on automatically

(The unit starts with "All in one", "Warmth: Low", "Intensity: Low" of the "Toe ~ Knee")

- The LED lamps of "All in one", "Warmth: Low" and "Intensity: Low" light after all the LED lamps light in turn.
- The warmer (heater) begins to warm after 5 minutes.  
(The warmer (heater) cannot operate independently without massaging)
- The power will be cut off after about 10 minutes.  
(The power will be cut off after about 10 minutes from the start, even when the program is switched over halfway)
- If you want to stop the massage halfway, press  repetitively until it turns to "Off"

## To change the program

- By pressing the key, the Program indicator (LED) will switch in the following order

Toe ~ Knee

→ All in one → Squeeze → Knead → Toe ~ Calf → Rear Knee

Toe ~ Thigh

→ All in one → Soft Knead → Deep Tissue


## Program content

Toe ~ Knee					Toe ~ Thigh		
All in one	Squeeze	Knead*	Toe ~ Calf	Rear Knee	All in one	Soft Knead*	Deep Tissue*
Specially selected massage with various techniques such as "Squeeze" and "Knead".	Squeeze your toes to the back of your knees like squeezing. Strong massage.	Knead the whole legs fully from the top to the bottom in turn.	Complex massage with "Squeeze" and "Knead" focusing on the toes to calves.	Complex massage with "Squeeze" and "Knead" focusing on the back of the knees.	Specially selected massage with various techniques such as "Squeeze" and "Knead".	Compress the whole of both legs softly like you are wrapping them up and knead them slowly and comfortably.	Knead and loosen the muscles on the front, back and outside of your thighs a little stronger.

\*Pressure function like knead

## WARNING

- For initial use, select "Intensity: Low" mode until you are accustomed to it.  
(Otherwise it may lead to injuries or poor body conditions.)

## Toe ~ Knee

Try it with your legs stretched if you feel less pressure of the point balls.

## To change the intensity


- By pressing the key, the Intensity indicator (LED) will switch in the following order

Low → Middle → High → Off

## Adjust temperature around toes

- By pressing the key, the Warmth indicator (LED) will switch in the following order

Low → High → Off

## Warmer (heater)

- The warmer (heater) warms the textile around toes only.

Massage

(To change the program) (Adjust temperature around toes) (To change the intensity)

## After use

Be sure that the power supply is in the "Off" state.

### Remove the attachments from your legs and unplug the AC adaptor / power plug / air tube plug


Operate according to the steps opposite to those on Page EN14-EN15, EN16-EN17

#### Storage

- 1 Fasten the surface adhesive tape on the toe, ankle and calf, and fold the toe inside
- 2 Fold the straps of ankle and calf, roll up the air tube for connection and fix with the surface adhesive tape of the strap

#### CAUTION

- When pulling out the air tube plug, do not pull the air tube for connection.
  - For storage, do not wrap the power cord around the manipulator or the power plug.
  - Do not fold the air tube for connection or wrap it on the attachment.
- (Otherwise, the power cord or air tube for connection may be disconnected due to an overload of the power cord, causing fire or electric shock due to short circuit.)

#### Keeping

- Do not place it within the reach of children
- Do not place it at a location with high temperature, high humidity or direct sunlight
- Do not place any objects on the manipulator and attachment surface  
(Otherwise the traces may remain on the surface or lead to cross-color or discoloration.)

## Cleaning and maintenance


Be sure that the power supply is in the "Off" state.

Do not clean the unit by thinner, benzine or alcohol.  
(Otherwise it may lead to malfunctioning, component damage or discoloration.)

- 1 Dip the cloth in water and wring out the water  
(For serious dirt, dilute a neutral cleaning agent with warm water or a water solution of about 5%.)

- 2 Wipe the dirt

- Please avoid using excessive force
- Do not immerse manipulator or attachment in water or wash with water.  
(Doing so may cause electric shock, burning or fire due to a short circuit.)

If the surface adhesive tape is stained with dirt or hair, clean it with a toothpick or other tools

- Do not put foreign objects into the insertion opening for air tube plug

# Frequently asked questions

## Under what circumstances can it not be used?

- Pregnant women and women who have just given birth
  - Persons suffering from severe circulatory disorders in the legs
- You cannot use the unit in the above cases. (Page EN4)
- If you have a fever or symptom in your body, consult doctor before use.

## Can I use it many times a day?

- It can be used once in a day. (The power will be cut off after about 10 minutes)
- If the time for a single use is less than 10 minutes, the total time for multiple uses shall be less than 10 minutes per day.
- After using it for a long time, it will produce excessive stimulation against the muscles and nerves and sometimes it will lead to the opposite effect.

### **WARNING**

- **The service time shall be less than 10 minutes per day**  
(Otherwise it may lead to a low-temperature burn, the opposite effect or injuries.)

## Can I use the unit for parts other than feet / calves / around knees / thighs?

- The unit is for feet / calves / around knees and thighs.
- If it is used for other parts, it may lead to accidents or injuries.

## What is the size of the foot / calf / thigh?

- Available foot size: approx. 20 cm ~ 28 cm
- Available calf girth: approx. 27 cm ~ 48 cm
- Available thigh girth: approx. 45 cm ~ 64 cm

## The toes are too hot...

- Set the warmth to "Low". If it is still hot, put on thick socks.

## The ankles and calves are not warm...

- The warmer (heater) warms the textile around toes only.

## The massage intensity is too high...

## The massage intensity is too low...

- The massage intensity can be adjusted by the wrapping method of the attachment.  
Too high intensity → wrap the attachment loosely.  
Too low intensity → wrap the attachment tightly.
- Try to use the sole and shin pads when the massage intensity for soles and shins is weak. (Page EN18-EN19)

## I feel uncomfortable when the skin is in contact with the surface adhesive tape...

- Cover the skin with a thin handkerchief or other article.

## The air in the attachment cannot be discharged completely...

- If there is a lot of air remaining in the attachment (airbag), it will take some time to discharge. Remove the air tube plug and wait a moment. The air can be discharged on its own.

## It sends out the following noise. Is it abnormal?

- Crackling sound  
(Tensioning sound of surface adhesive tape when the air enters into the attachment (airbag).)
  - Kachi sound  
(It is the sound sent out when the manipulator adjusts the air.)
  - Humming sound  
(Motor sound.)
- It is the sound sent by the unit structure. This is normal.

## Can it be used abroad?

The unit adopts the power supply of 220 V – 240 V ~.

In addition, when using a transformer or a conversion plug, a malfunction may occur. Do not use them.

# Before entrusted repair

## Please confirm the following.

After confirmation, if the abnormality remains, first unplug the AC adaptor from the socket and then entrust to an authorized service center for inspection and repairs.

Symptoms	Items to be confirmed
The power cannot be switched on	<ul style="list-style-type: none"> <li>Does the AC adaptor fall off? (Page EN14, 16)</li> </ul>
It stops halfway	<ul style="list-style-type: none"> <li>Whether the Power / Intensity selection key is in the "Off" state? (Page EN20-EN21)</li> <li>Does the air tube plug fall off? (Page EN15, 17)</li> </ul>
The air cannot enter	<ul style="list-style-type: none"> <li>Whether the air tube for connection is folded?</li> </ul>
Toes do not feel warmth	<ul style="list-style-type: none"> <li>The warmer (heater) warms in 5 minutes after turning on the unit.</li> <li>Is the Warmth indicator (LED) on? Pressing the Warmth key lights the Warmth indicator (LED). (Page EN21)</li> <li>Are your heels adjusted to the specified position? (Page EN15, 17)</li> <li>Are you using the sole pads? (Page EN18)</li> <li>Are you putting the unit on thick socks?</li> </ul>
I feel warmth differently on right and left	<ul style="list-style-type: none"> <li>Does one air tube plug fall off? (Page EN15, 17)</li> <li>Are your heels adjusted to the specified position? (Page EN15, 17)</li> </ul>

Symptoms	Items to be confirmed
The massage intensity is low	<ul style="list-style-type: none"> <li>Is the attachment wrapped loosely? (Page EN15, 17)</li> <li>Are you putting the unit on thick socks?</li> <li>Does the air tube plug fall off? (Page EN15, 17)</li> <li>Is the air tube plug inserted completely? (Page EN15, 17)</li> <li>Whether the air tube for connection is folded?</li> </ul>
The massage intensity is high	<ul style="list-style-type: none"> <li>Is the attachment wrapped tightly? (Page EN15, 17)</li> </ul>
Intensity indicator (LED), Program indicator (LED) or Warmth indicator (LED) does not go on	<ul style="list-style-type: none"> <li>Does the AC adaptor fall off? (Page EN14, 16)</li> <li>Press the Power / Intensity selection key, Program selection key and Warmth key and confirm. (Page EN20-EN21)</li> </ul>
The air cannot be discharged completely	<ul style="list-style-type: none"> <li>Remove the attachment and unplug the air tube plug. (Page EN22)</li> </ul>

## In case of abnormalities ...

Manipulator or AC adaptor is overheating	<p>These conditions are failures. To prevent an accident, stop using immediately and contact the dealer.</p>
It is damaged	
The attachment cloth is damaged	

Before entrusted repair

# Specifications

<b>Product name</b>		Air Massager
<b>Product model</b>		EW-RAJ6
<b>Power supply</b>		220 V – 240 V ~ 50 – 60 Hz
<b>Power consumption</b>		24 W (Rated power of the heat transfer product 5 W)
<b>Air pressure</b>		Low : approx. 36 kPa Middle : approx. 44 kPa High : approx. 54 kPa (It slightly varies depending on the program or part)
<b>Duration</b>		Approx. 10 minutes
<b>Dimension (packing qty)</b>	<b>Manipulator (1 pc)</b>	Approx. 145 mm (D) × 60 mm (H)
	<b>Attachment (2 pcs)</b>	Approx. 780 mm (W) × 300 mm (L) × 810 mm (H)
<b>Weight</b>		Approx. 2 390 g (manipulator + attachment)
<b>Length of power cord</b>		Approx. 1.8 m
<b>Corresponding perimeter of attachment</b>		Foot size : approx. 20 cm ~ 28 cm Calf girth : approx. 27 cm ~ 48 cm Thigh girth: approx. 45 cm ~ 64 cm
<b>Surface temperature of the inside textile of the attachment (in indoor temperature of 20 °C)</b>		Low: approx. 37 °C High: approx. 40 °C
<b>Manufacturer</b>		Panasonic Manufacturing (Shanghai) Co., Ltd. No. 258 East Jiangtian Road, Songjiang Industrial Park, Shanghai
<b>Country of origin</b>		China

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