

Bathroom Equipment



Product Details


Useful products for maintenance of bathroom environment!

Yukuria series with antimicrobial processing

Antimicrobial mechanism

Making an antimicrobial agent and antimicrobial barrier

Blocking the growth of mildew inside



※ Antimicrobial processing cannot eliminate mildew. Mildew may occur depending on conditions of use.

Can easily rotate when sitting on shower chair

Rotatory Shower Chair

Maximum user weight 100Kg



Orange

Safe, comfortable and soft cushion



The underside of the armrests is soft too. No worry even when the user accidentally bumps into the armrests when transferring to the bathtub.



30CM

- The soft seat surface also let users of skinny buttocks sit comfortably.
- ※ This material does not let you easily feel the coolness and is not slippery regardless of being in wet or dry conditions.
- Easy to keep dry and removal
- ※ If soap or shampoo attached on the surface, it would become slippery and may cause injury due to falling.

It supports the body when getting in or out of the bathtub

Bath Grip

Assembled finished product
Leg attached type only.
Leg mounting is required.

Maximum user weight 100Kg

Applicable for a wide range of bathtubs
Bath Grip (specially for bathtub constructed by conventional method)

For UB only
※ with compact Bath Grip leg (PN-L12111D)



For unit bath, do not push in the apron part.
For UB only ※ bath grip!

As it is applicable for a wide range of bathtub curves, it can be pushed to the edge.



If it is pushed to the edge, it is easy to step over the bathtub.



Please check the details in advance for tub width, curves etc.

To make up the height difference of the bathtub

Small Bath Stool


Maximum user weight 100Kg

Small Bath Stool



Easy to set up

As it sinks under its own weight, it is easy to set up even after the bathtub has been filled with hot water.



Mini suckers which are slip resistant even in the bathroom and inside the bathtub have been adopted.

Even there is a change in the way that the user gets on his feet, the height is adjustable.



In order to prevent dislocation after the surgery of artificial hip joint due to femoral neck fracture, the height of the seat surface can be set higher.

Soft type is more comfortable for sitting
The detachable cover can be washed



Soft touch to foot
Soft touch to thighs



* Please consult your doctor or caregiver if necessary.